M40 Colonoscope Training Model

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

Case 6. ("Reverse alpha" loop formation) p.16 • p.17

Preparation and Set-up


First, attach the fixtures to the colon-rectum tube in order, as illustrated bellow, threading the tube through the fixture rings. Keep in mind the two types of fixtures, Type A in positions 2, 3 and 8, and Type B for the rest.


After filling the colon-rectum tube with lubricant (see pages 3 & 5 of the "Instruction Manual"), attach it to the anus opening of the abdomen unit. Close the opening with the supplied plug after affixing the plate.



- Place one of the colon layout guide boards at the bottom of the abdomen model.

- Take off the four screws from the white plate and remove the plate. Insert the anal end of the colon-rectum tube through the opening of the abdomen model wall.

- Fit the flange at the tube-end to the black base without pinching or creasing.

- Attach the white plate to the base so that the manual pump air tube fits in the notch. Affix the plate with the four screws.


Primary level

Case 1

■ Introductory level I

■ Simple layout appropriate for first-time trainees.

■ Objectives


Learn how to insert the colonoscope deeply into the transverse colon and the ascending colon, without forming a loop at the sigmoid colon.

Acquire basic insertion skills required to pass through at each part of the colon.

■ Skills Gained

Learn the "hooking the fold" method to pass through the sigmoid colon.

Learn "withdrawal" manipulation to go through the hepatic flexure.


Primary level

Case 2

■ Introductory level II

- Still a simple layout, but with a longer sigmoid colon than in Case 1.
An essential training step before moving on to actual patients.


■ Objectives

Pass through the sigmoid-descending colon (S-D) junction without forming a loop, while practicing a right withdrawal with clockwise torque.

■ Skills Gained

Advance the colonoscope into the descending colon without forming a loop while practicing the "hooking the fold" method with sufficient suction and straightening of the scope.


Primary level

Case 3

■ Straighten out the "Alpha" loop

- Still a simple layout, but with a naturally formed "Alpha" loop in a longer sigmoid colon. With insertion, the colonoscope naturally forms an "alpha" loop, and advances toward the splenic flexure.


■ Objectives

Learn how to straighten out an "alpha" loop at the sigmoid colon.

■ Skills Gained

Reduce the "Alpha" loop at the sigmoid colon, practicing the "right turn shortening" method (withdrawal with clockwise torque).


Secondary level

Case 4

■ Redundant sigmoid colon, "N" loop formation and long transverse colon

■ A challenging layout, with a difficult-to-shorten sigmoid colon and a drooping transverse colon.


■ Objectives

Manage successful insertion for challenging cases with a long sigmoid colon and drooping transverse colon.

■ Skills Gained

Proper insertion techniques utilizing all skills and methods learned at the primary level.

To accomplish even higher skill, trainees may be able to challenge to insert the colonoscope without forming any loop, which is considered to give less pain to patients.


Secondary level

Case 5

■ Redundant sigmoid colon and "alpha" loop formation

- An "alpha" loop is formed by the redundant sigmoid colon, therefore, passing through the sigmoid colon by shortening method is extremely challenging.


■Objectives

Learn how to reduce a sigmoid colon loop and how to achieve the "shortening of the colon" when trainee encounter challenging cases.

■Skills Gained

Proper insertion techniques utilizing all skills and methods learned at the primary level.

To accomplish even higher skill, trainees may be able to challenge to insert the colonoscope without forming any loop, which is considered to give less pain to patients.

Case 6

- Because of the significantly redundant sigmoid colon, a "reverse alpha" loop is formed by insertion of the colonoscope.


■ Objectives

Acquire skills in reduction of a "reverse alpha" loop, typically seen in lean patients, by using left torque.

■ Skills Gained

Since the shortening of the sigmoid colon is difficult to achieve, insert and advance the instrument first, until the tip reaches to the descending colon. Then straighten out the "reverse alpha" loop by left torque without letting the tip of the colonoscope drop off when it is pulled back.

KYOTO KAGAKU co.,LTD

URL: <http://www.kyotokagaku.com> e-mail: rw-kyoto@kyotokagaku.co.jp

Worldwide Inquiries & Ordering

Kyotokagaku Head Office and Factories:

15 Kitanechoya-cho, Fushimi-ku, Kyoto, 612-8388, JAPAN

Tel: +81-75-605-2510 Fax: +81-75-605-2519

All American regions

KYOTOKAGAKU AMERICA INC.

3109 Lomita Boulevard, Torrance, CA 90505-5108, USA

Tel: 1-310-325-8860 Fax: 1-310-325-8867

The contents of the instruction manual are subject to change without prior notice.

No part of this instruction manual may be reproduced or transmitted in any form without permission from the manufacturer. Please contact manufacturer for extra copies of this manual which may contain important updates and revisions.

Please contact manufacturer with any discrepancies, typos, or mistakes in this manual or product feedback. Your cooperation is greatly appreciated.

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